# HEMATOLOGY CONTROL

## MEK-3DL, MEK-3DN, MEK-3DH

English

### General

#### Intended Purpose

*For in vitro diagnostic use only.*

HEMATOLOGY CONTROL is for in vitro diagnostic use as a value assigned whole blood CONTROL designed to monitor the values of Nihon Kohden hematology analyzers. HEMATOLOGY CONTROL is not used for diagnosis or assistance in diagnosis, prognosis, prediction of disease or patient's physiological condition. Measuring this HEMATOLOGY CONTROL ensures that the device is functioning according to its intended use.

Control parameters are below:

- 1) WBC: White Blood Cell Count
- 2) RBC: Red Blood Cell Count
- 3) HGB: Hemoglobin Concentration
- 4) HCT: Hematocrit Percent
- 5) MCV: Mean Corpuscular Volume
- 6) MCH: Mean Corpuscular Hemoglobin
- 7) MCHC: Mean Corpuscular Hemoglobin Concentration
- 8) PLT: Platelet Count
- 9) MPV: Mean Platelet Volume
- 10) RDW: Red Blood Cell Distribution Width
- 11) LY%: Lymphocyte Percent
- 12) MO%: Monocyte Percent
- 13) GR%: Granulocyte Percent
- 14) LY: Lymphocyte
- 15) MO: Monocyte
- 16) GR: Granulocyte
- 17) QC ESR: QC Erythrocyte Sedimentation Rate

Read the analyzer operator's manual together with this manual before and during use.

Also refer to the assay sheet provided with HEMATOLOGY CONTROL.

**NOTE** • Use HEMATOLOGY CONTROL with the specified analyzers only.

- Read the SDS (Safety Data Sheet) carefully before use. The SDS is available from your Nihon Kohden representative.

#### Devices Intended for Use in Combination with

- MEK-1301 automated hematology analyzer
- MEK-1302 automated hematology analyzer
- MEK-1305 automated hematology and ESR analyzer
- MEK-6400 automated hematology analyzer
- MEK-6410 automated hematology analyzer
- MEK-6420 automated hematology analyzer
- MEK-6450 automated hematology analyzer (for animal)
- MEK-6500 automated hematology analyzer
- MEK-6510 automated hematology analyzer
- MEK-6550 automated hematology analyzer (for animal)

#### Materials Provided and Materials Required

##### Materials Provided

HEMATOLOGY CONTROL

##### Materials Required (Not Provided with HEMATOLOGY CONTROL)

- Protective gloves
- Protective glasses
- Waste container for biohazardous materials

#### Intended Users

For laboratory professional use only, in laboratories with suitable equipment for hematological testing. Qualified personnel, e.g. laboratory technicians trained in hematology analysis techniques, will be able to use according to the operator's manual.

#### Symbols

The following symbols are used with HEMATOLOGY CONTROL. The descriptions of each symbol are given in the table below.

Symbol	Description	Symbol	Description
	Caution		Operator's manual; operating instructions

0614-908006

Symbol	Description
	Catalogue number
	Lot number
	Use by
	In vitro diagnostic medical device
	Manufacturer
	European representative
	Temperature limits
	Biological risks
	This way up

Symbol	Description
	Control/ Concentration level
	Open vial stability
	Unique Device Identifier
	The CE mark is a protected conformity mark of the European Union. The four digits after the CE mark indicate the identification number of the Notified Body involved in assessing the product's conformity as a medical device.

## Safety Information

**⚠ WARNING** A warning alerts the user to the possible injury or death associated with the use or misuse of the instrument.

**⚠ CAUTION** A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

Pay attention to all safety information in this operator's manual.

### ⚠ WARNING

- **POTENTIALLY BIOHAZARDOUS MATERIAL.** For in vitro diagnostic use. Each human donor/unit used in the preparation of HEMATOLOGY CONTROL has been tested by an FDA licensed method/test and found to be negative or non-reactive for the presence of HBsAg, Anti-HCV, NAT testing for HIV-1, HCV (RNA) and HIV-1/2. Each unit is also negative by a serological test for Syphilis (RPR or STS). Because no test method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious. When handling or disposing of vials follow precautions for patient specimens as specified in the OSHA Bloodborne Pathogen Rule (29 CFR Part 1910, 1030) or other equivalent biosafety procedures.
- Wear protective equipment such as disposable gloves when handling HEMATOLOGY CONTROL.
- Do not swallow HEMATOLOGY CONTROL.

### ⚠ CAUTION

HEMATOLOGY CONTROL is intended for in vitro diagnostic use only by trained personnel.

- NOTE**
- Do not use HEMATOLOGY CONTROL if it is past the expiration date on HEMATOLOGY CONTROL or if HEMATOLOGY CONTROL is stored under unspecified conditions.
  - Store HEMATOLOGY CONTROL in a pharmaceutical refrigerator.
  - Store HEMATOLOGY CONTROL at temperatures between 2 and 8°C (36 and 46°F).
  - Do not freeze HEMATOLOGY CONTROL.
  - If the cap of HEMATOLOGY CONTROL becomes dirty, wipe with a non-woven towel such as Kimwipes.
  - Platelet cannot be measured by precipitation or centrifuge.
  - This product does not support measuring WBC differential under the microscope.

## Using HEMATOLOGY CONTROL

- NOTE**
- Refer to the analyzer operator's manual for details on measurement.
  - Do not use HEMATOLOGY CONTROL if deterioration is suspected. After mixing, HEMATOLOGY CONTROL should be similar in appearance to fresh whole blood. In unmixed tubes, the supernatant may appear cloudy and reddish; this is normal and does not indicate deterioration. Other discoloration, very dark red supernatant or unacceptable results may indicate deterioration.
  - Confirm that the lot number of the product matches the lot number on the assay sheet.
  - When using the MEK-1301 or MEK-1302 automated hematology analyzer or MEK-1305 automated hematology and ESR analyzer, read the three barcodes printed on the "Barcode Entry" area on the assay sheet provided with HEMATOLOGY CONTROL with the connected barcode reader. The barcodes may be read in any order.

Measurement of a control is recommended for statistical quality control of the measurement accuracy of the analyzer. Also, it is recommended to use the average and upper and lower limits set by each laboratory to measure the control.

## Frequency of Performing Quality Control

Determine according to the procedures established by each laboratory.

## Control Used

HEMATOLOGY CONTROL

## Measurement Principle

Refer to the analyzer operator's manual.

## Procedure

1. Take HEMATOLOGY CONTROL out of the pharmaceutical refrigerator and allow it to warm to room temperature (15 to 30°C, 59 to 86°F) for 15 minutes before mixing.

2. Mix HEMATOLOGY CONTROL.


**NOTE**

- Do not use a mechanical mixer.
- Do not mix HEMATOLOGY CONTROL roughly. If it bubbles, hemolysis is caused.

- 1) Hold a tube horizontally between the palms of the hands and roll it back and forth for 20 to 30 seconds. Occasionally invert the tube so that HEMATOLOGY CONTROL is mixed thoroughly.

Roll the tube back and forth

Invert the tube

- 2) Continue to mix in the manner described in step 1) until the red cells are completely suspended. Tubes stored for a long time may require extra mixing.
3. Gently invert the tube until there is no clot of blood on the bottom before measurement.
4. Measure HEMATOLOGY CONTROL as instructed in the analyzer operator's manual.
5. After measurement, clean residual material from the cap and tube rim with a lint-free tissue if the tube has been open for measurement. Replace the cap tightly.
6. Return the tube to the refrigerator within 30 minutes.

## Technical Information

### Metrological Traceability of Values Assigned to HEMATOLOGY CONTROL

For each parameter, the internationally recognized reference method is:

- 1) WBC and RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT and PCT <sup>1</sup>: CLSI H7-A3
- 4) PLT: ICSH/WHO 2000: Recommended methods for the visual determination of white cell count and platelets count. WHO/DIL/00.3, 2000

5) LY, MO and GR: CLSI H26-A2:Vol.29 No.40 and CLSI H20-A2:Vol.27 No.4

6) ESR: CLSI H02-A5

<sup>1</sup> PCT: Platelet Crit

## Expected Values

- Verify that the lot number on the tube matches the lot number on the assay sheet.
- Refer to the assay values for the specific instrument model on the provided assay sheet.
- Assay values are determined on well-maintained, properly calibrated instruments using the conformed reagents. Reagent differences, maintenance, operation technique, and calibration may contribute to inter-laboratory variation.

**NOTE:** RDW-CV is displayed as RDW on the MEK-6400, MEK-6410, MEK-6420 and MEK-6450 automatic hematology analyzers.

## Composition

HEMATOLOGY CONTROL contains human erythrocytes, simulated leukocytes and mammalian platelets suspended in a plasma-like fluid with preservatives.

## Sterilization Method

HEMATOLOGY CONTROL is not intended to be sterilized or kept in a sterile environment.

## Interfering Substances or Limitations

The performance of the control is assured only if it is properly stored and used as described in this operator's manual. Incomplete mixing of the tube prior to use invalidates both the sample withdrawn and any remaining material in the tube.

## Analytical Performance Characteristics

- Assigned values are presented as a Mean and Range. The Mean is derived from replicate testing on instruments operated and maintained according to the manufacturer's instructions. The Range is an estimate of variation between laboratories and also takes into account inherent imprecision of the method and expected biological variability of the control material.
- Assay values on a new lot of control should be confirmed before the new lot is put into routine use. Test the new lot when the instrument in good working order and quality control results on the old lot are acceptable. The laboratory's recovered mean should be within the assay range.
- For greater control sensitivity each laboratory should establish its own mean and acceptable range and periodically reevaluate the mean. The laboratory range may include values outside of the assay range. The user may establish assay values not listed on the assay sheet, if the control is suitable for the method.

## Mathematical Approach upon which the Calculation of the Analytical Result is Made

The following parameters are calculated as follows <sup>1</sup>:

$$MCV = HCT / RBC$$

$$MCH = HGB / RBC$$

$$MCHC = HGB / HCT$$

$$MPV = PCT / PLT$$

$$LY\% = LY / WBC$$

$$MO\% = MO / WBC$$

$$GR\% = GR / WBC$$

<sup>1</sup> Without unit correction

- NOTE
- RBC indices are calculated by rounding off to the whole number or the nearest tenth.
  - QC ESR is not an actual erythrocyte sedimentation rate. The value was calculated specifically for quality control of a Nihon Kohden analyzer based on the optical density, HCT and MCV of HEMATOLOGY CONTROL.

## Environmental Conditions

### Storage and Transport Environment

Temperature: 2 to 8°C (36 to 46°F)

### Usage Environment

Temperature: 15 to 30°C (59 to 86°F)

## Expiration Date

The expiration date is shown on the container label.

Lifetime: 3 months

## Shelf Life Date After Opening

14 days (within the expiration date)

## Package and Catalog Number

Model	Qty	Catalog Number
MEK-3DN	MEK-3DN (2.0 to 2.4 mL) × 3	3DN03
	MEK-3DN (2.0 to 2.4 mL) × 6	3DN06
	MEK-3DN (2.0 to 2.4 mL) × 12	3DN12
MEK-3DL, MEK-3DH	MEK-3DL (2.0 to 2.4 mL) × 3, MEK-3DH (2.0 to 2.4 mL) × 3	3DLH3
	MEK-3DL (2.0 to 2.4 mL) × 6, MEK-3DH (2.0 to 2.4 mL) × 6	3DLH6

## Disposal

### ⚠ WARNING

Dispose of HEMATOLOGY CONTROL according to your local laws and your facility's guidelines (including incineration, melt treatment, sterilization, disinfection and request for waste disposal) for disposing of infectious medical waste. Otherwise, it may affect the environment. If there is a possibility that HEMATOLOGY CONTROL may have been contaminated with infection, it may cause infection.

When disposing of HEMATOLOGY CONTROL, such as when the expiration date is past, follow the instructions on the SDS of HEMATOLOGY CONTROL.

## Revision History

Edition	Date	Details	Code Number
1st Edition	07 Mar 2016	Initial issue	0604-910978
4th Edition	30 Mar 2023	IVDR compliance	0614-908006

- NOTE
- The code number of this manual was changed from 0604-910978B to 0614-908006 when the manual was updated from 3rd Edition to 4th Edition.
  - Changes made in the most recent edition are indicated by a bar in the left margin of each page.

Note for users in the territory of the EEA and Switzerland:  
Any serious incident that has occurred in relation to the device must be reported to the European Representative designated by the manufacturer and the Competent Authority of the Member State of the EEA and Switzerland in which the user and/or patient is established.

### Copyright Notice

The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved.

 Manufacturer  
**NIHON KOHDEN CORPORATION**  
1-31-4 Nishiochiai, Shinjuku-ku,  
Tokyo 161-8560, Japan  
Phone +81 3-5996-8041  
<https://www.nihonkohden.com/>

 European Representative  
**NIHON KOHDEN EUROPE GmbH**  
Raiffeisenstrasse 10, 61191 Rosbach, Germany  
Phone +49 6003-827-0 Fax +49 6003-827-599  
UK Responsible Person  
**NIHON KOHDEN UK LTD.**  
Unit 3, Heyworth Business Park,  
Old Portsmouth Road, Peasmarsh,  
Guildford, Surrey, GU3 1AF, UK  
Phone +44 14 8333-1328

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