

Transmitters

ZS/ZM series

Better care by **wireless**



Fighting Disease with Electronics

 **NIHON KOHDEN**

High quality patient care by wireless

— patient is free from the bed —

Wireless monitoring system

You can make a wireless monitoring system with transmitters and a Nihon Kohden central monitor with multiple patient receiver or a Nihon Kohden telemetry system.

Increase patient safety

By continuous measurement of important vital parameters, patient condition is always monitored so changes can be quickly handled.

Increased patient comfort

No cable between the patient and monitor, so the patient is free from the bed. The transmitter is compact and light-weight, so it is easy to wear.

More efficient patient care

Instead of visit and manual measurement, automatic measurement saves the care staff's time.

Easy to use

No complicated preparation to start measurement. The CALL key activates the recorder* and nurse call at the receiving monitor (ZS-611/620/630/530P). The long battery life reduces running costs.

*Depends on setting


Related Products

Central monitor
CNS-6201

Central monitor
CNS-9101

Telemetry system
WEP-5204/5208

Multiple patient receiver
ORG-9100



Measurement Parameters



ZS-611P

SpO₂
5.5 days*
2 AA (LR6)
batteries



ZS-620P

ECG, Resp
6 days*
1 AA (LR6)
battery



ZM-920P

ECG(8-lead),
Resp
8 days*
3 AA (LR6)
batteries



ZS-630P

ECG, Resp, SpO₂
4 days*
2 AA (LR6)
batteries



ZS-530P

ECG, Resp, SpO₂
3 days*
2 AA (LR6) batteries
• 2.2-inch TFT color LCD
• 10 minutes review data



ZM-930P

ECG(8-lead),
Resp, SpO₂
5 days*
2 AA (LR6)
batteries



ZS-940P/ZM-940P

NIBP, ECG, Resp, SpO₂
2 days*
3 AA (R6) NiMH batteries
• LCD shows NIBP, SpO₂ and pulse rate
• NIBP measurement mode can be selected from manual, continuous and interval
• 8-lead ECG measurement is available (ZM-940P)

* under certain conditions

This brochure may be revised or replaced by Nihon Kohden at any time without notice.

NIHON KOHDEN

NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku, Tokyo 161-8560, Japan
Phone +81 (3) 5996-8036 Fax +81 (3) 5996-8100
www.nihonkohden.com