

# External Paddle Holder

## DP-831VK

### General

The DP-831VK external paddle holder is installed on a TEC-8321K, TEC-8322K, TEC-8332K, TEC-8342K and TEC-8352K defibrillator and lets you put the optional ND-831V external paddles on it.

### Safety Information

**CAUTION** A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

Pay attention to all safety information in the Operator's Manual or Installation Guide.

#### CAUTION

Before installing or removing the holder, make sure that the defibrillator is turned off and the AC/DC module and the battery pack are removed. Otherwise, the operator may receive electrical shock or injury.

This Safety and Performance Information is an extract from the general and safety information sections of the most recent edition of Operator's Manual or Installation Guide. Therefore, the contents of your Operator's Manual or Installation Guide may differ from those of this Safety and Performance Information. For detailed operating procedures, follow the instructions of your Operator's Manual or Installation Guide.

 Manufacturer  
**NIHON KOHDEN CORPORATION**  
1-31-4 Nishiochiai, Shinjuku-ku,  
Tokyo 161-8560, Japan  
Phone +81 3-5996-8041  
<https://www.nihonkohden.com/>

 European Representative  
**NIHON KOHDEN EUROPE GmbH**  
Raiffeisenstrasse 10, D-61191 Rosbach, Germany  
Phone +49 6003-827-0 Fax +49 6003-827-599