

AUXOUT Unit

QI-831V

General

The QI-831V AUXOUT Unit connects a bedside monitor to a TEC-8321K, TEC-8322K, TEC-8332K, TEC-8342K or TEC-8352K defibrillator and transfers ECG measured with the defibrillator to the bedside monitor.

Safety Information

WARNING A warning alerts the user to possible injury or death associated with the use or misuse of the instrument.

CAUTION A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

Pay attention to all safety information in the Operator's Manual or Installation Guide.

WARNING

Connect the AUXOUT unit to the specified defibrillator and follow the specified procedure. Failure to follow this warning may result in electrical shock or injury to the patient and operator.

CAUTION

Before installing or removing the unit, make sure that the defibrillator is turned off and the AC/DC module and the battery pack are removed. Otherwise, the operator may receive electrical shock or injury.

This Safety and Performance Information is an extract from the general and safety information sections of the most recent edition of Operator's Manual or Installation Guide. Therefore, the contents of your Operator's Manual or Installation Guide may differ from those of this Safety and Performance Information. For detailed operating procedures, follow the instructions of your Operator's Manual or Installation Guide.

 Manufacturer
NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>

 European Representative
NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, D-61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599