


MEK-5DL HEMATOLOGY CONTROL
for WBC 5 part differential 3.0mL

ASSAY SHEET

LOT	MEK5D23YL
	2024-02-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±0.6	3.0	3.0	3.0	3.00
NE%	%	±15.0	28.7	28.7	28.7	28.46
LY%	%	±15.0	53.9	53.9	53.9	57.58
MO%	%	±5.0	11.9	11.9	11.9	9.43
EO%	%	±5.0	2.5	2.5	2.5	2.96
BA%	%	±4.0	3.0	3.0	3.0	1.57
NE	10 ³ /μL	±0.5	0.8	0.8	0.8	0.85
LY	10 ³ /μL	±0.5	1.6	1.6	1.6	1.73
MO	10 ³ /μL	±0.2	0.4	0.4	0.4	0.28
EO	10 ³ /μL	±0.2	0.1	0.1	0.1	0.09
BA	10 ³ /μL	±0.2	0.1	0.1	0.1	0.05
RBC	10 ⁶ /μL	±0.15	2.62	2.61	2.63	2.58
HGB	g/dL	±0.4	6.5	6.6	6.5	6.62
HCT	%	±2.0	21.7	21.0	20.6	23.0
MCV	fL	±7.0	82.9	80.3	78.2	89.1
MCH	pg	±2.4	24.8	25.3	24.7	25.7
MCHC	g/dL	±3.0	30.0	31.4	31.6	28.8
RDW-CV *	%	±4.0	15.8	17.5	16.2	16.1
PLT	10 ³ /μL	±30	57	64	57	60.0
MPV	fL	±3.0	4.4	5.4	4.7	7.2

Manual Registration	
Lot	23Y
Level	Low
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901690M