


MEK-5DH HEMATOLOGY CONTROL
for WBC 5 part differential 3.0mL

ASSAY SHEET

LOT	MEK5D242H
	2024-05-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±2.2	24.9	24.8	24.4	25.86
NE%	%	±10.0	67.9	67.9	67.9	68.02
LY%	%	±10.0	17.2	17.2	17.2	18.24
MO%	%	±5.0	9.4	9.4	9.4	7.20
EO%	%	±5.0	4.5	4.5	4.5	5.65
BA%	%	±4.0	1.0	1.0	1.0	0.89
NE	10 ³ /μL	±2.8	17.0	16.9	16.6	17.59
LY	10 ³ /μL	±2.8	4.3	4.3	4.2	4.72
MO	10 ³ /μL	±1.0	2.3	2.3	2.3	1.86
EO	10 ³ /μL	±1.0	1.1	1.1	1.1	1.46
BA	10 ³ /μL	±0.8	0.2	0.2	0.2	0.23
RBC	10 ⁶ /μL	±0.25	5.49	5.45	5.44	5.64
HGB	g/dL	±0.7	17.3	17.1	17.2	17.36
HCT	%	±4.0	53.7	52.2	50.7	61.3
MCV	fL	±8.0	97.8	95.7	93.2	108.6
MCH	pg	±2.8	31.5	31.4	31.6	30.8
MCHC	g/dL	±3.0	32.2	32.8	33.9	28.3
RDW-CV *	%	±4.0	13.5	14.3	13.7	14.9
PLT	10 ³ /μL	±70	598	588	596	596.6
MPV	fL	±4.0	8.1	8.3	8.3	10.2

Manual Registration	
Lot	242
Level	High
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901692N