



MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS244
	2024-05-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	9.1	9.1	9.3	9.2	9.0	9.30	9.08
RBC	10 ⁶ /μL	4.64	4.64	4.64	4.61	4.63	4.69	4.64
HGB	g/dL	14.0	14.1	13.8	13.8	13.7	14.03	13.98
HCT	%	40.8	40.8	42.5	41.3	40.1	47.4	46.8
MCV	fL	88.0	88.0	91.5	89.5	86.6	101.1	100.8
MCH	pg	30.2	30.4	29.7	29.9	29.6	29.9	30.1
MCHC	g/dL	34.3	34.6	32.5	33.4	34.2	29.6	29.9
PLT	10 ³ /μL	245	245	239	242	241	253.9	260.2
MPV	fL	6.4	8.6	6.4	7.4	6.4	9.1	9.1
RDW-CV*	%	15.8	15.8	14.7	16.0	15.0	15.2	15.5

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	169
FL	85
SD	106
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901780N