


MEK-3DN HEMATOLOGY CONTROL 2.0mL

ASSAY SHEET

LOT	B245N
	2024-08-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order. If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.8	7.7	7.7	7.71
RBC	10 ⁶ /μL	±0.20	4.71	4.71	4.72
HGB	g/dL	±0.5	14.0	14.1	14.01
HCT	%	±3.0	40.0	40.0	45.8
MCV	fL	±7.0	84.9	84.9	97.0
MCH	pg	±2.4	29.7	29.9	29.7
MCHC	g/dL	±3.0	35.0	35.3	30.6
PLT	10 ³ /μL	±50	260	260	279.2
MPV	fL	±2.0	5.6	8.1	8.3
RDW-CV *	%	±3.0	16.7	16.7	16.4
LY%	%	±10.0	39.8	39.8	36.13
MO%	%	±5.0	5.0	5.0	1.98
GR%	%	±10.0	55.2	55.2	61.89
LY	10 ³ /μL	±1.0	3.1	3.1	2.79
MO	10 ³ /μL	±0.4	0.4	0.4	0.15
GR	10 ³ /μL	±1.0	4.2	4.2	4.77

MEK-1305
QC ESR
10±5 mm

Manual Registration	
Lot	245
Level	Normal
Control	MEK-3D

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



CE
0123 **IVD**
0694-901685P