


MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS247
	2024-08-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry




		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	9.2	9.2	9.4	9.3	9.2	9.43	9.25
RBC	10 ⁶ /μL	4.65	4.65	4.60	4.57	4.59	4.64	4.61
HGB	g/dL	13.5	13.6	13.3	13.2	13.2	13.59	13.66
HCT	%	40.0	40.0	41.0	39.8	38.7	45.2	45.0
MCV	fL	86.0	86.0	89.2	87.1	84.4	97.5	97.6
MCH	pg	29.0	29.2	28.9	28.9	28.8	29.3	29.6
MCHC	g/dL	33.8	34.0	32.4	33.2	34.1	30.1	30.4
PLT	10 ³ /μL	256	256	246	250	249	261.8	271.6
MPV	fL	6.6	8.7	6.6	7.5	6.6	9.2	9.2
RDW-CV*	%	16.5	16.5	15.1	16.6	15.4	16.3	16.2

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	170
FL	84
SD	112
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.

 <https://www.nihonkohden.com/spi/index.html>

 **NIHON KOHDEN CORPORATION**
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


MEK-CAL


PLUS247


2024-08-05

 
0123
0694-901780Q